amiea

product data sheet

Most suitable for mobile use

The amiea Linelle II is lightweight and very easy to operate and at the same time providing all the essentials. This basic, reliable device for micropigmentation is perfect for professionals who offer mobile services.



APPLICATIONS

Cosmetic micropigmentation for:

- areola / body
- eyebrows, eye-liners, lip-liners and / or filling in lip colour

FEATURES	BENEFITS
Adjustable needle speed	 High power Perfect pain management Minimal downtime 50 - 140 pps
Needle length adjustable	 Appropriate needle length for each area of application Range of depth 0.00 - 0.098 inches / 0.0 - 2.5 mm
Precision and power	 Low epidermal trauma Low noise and vibration Superior results
Values	• German technology

• R&D and production in Berlin

amiea

PERFECT RESULTS

Micropigmentation can be used to restore confidence and improve your overall appearance. It can be applied to eyebrows, eyelines and lips.



eyebrows before





eyebrows after



lips before

lips after

TECHNICAL PARAMETERS

Technical parameters

•	Туре:	A-9.9-140
•	Rated voltage:	15 V - (DC)
•	Power input:	7 V A max
•	Working frequency:	Selectable from 50 to 140 hits per second
•	Drive:	Precision DC motor
•	Operating mode:	Continuous operation

Operating conditions

•	Ambient temperature:	+50° to +77°F / 15° - 25°C
•	Relative humidity:	30% to 75%
•	Width x height x depth:	3.94 x 1.77 x 5.51 inches / 10.0 x 4.5 x 14.0 cm
•	Weight of handpiece:	1.76 oz / 50 g
•	Total weight:	0.88 lbs / 14.11 oz / 400 g
•	Manufactured in:	Germany according to quality standard ISO 9001

Patents

- US 6,505,530; US 6,345,553; US 8,029,527; US 9,504,814 EP 1 495 782; EP 1 618 915
- Additional patents pending